

NEOGARD® Quick Notes

Pot Life & Recoat Time with Accelerators



Date: January 7, 2010
From: Technical Department

70410

Accelerator	Amount	Pot Life @40°F/60%RH	Pot Life @75°F/50%RH	Pot Life @95°F/50%RH	Recoat Time @40°F/60%RH	Recoat Time @75°F/50%RH	Recoat Time @95°F/50%RH
none	none	> 7 days	> 7 days	> 7 days	48 hrs	14 hrs	12 hrs
7931	4 oz/5 gal	> 24 hrs	> 24 hrs	> 24 hrs	22 hrs	10 hrs	8 hrs
7931	8 oz/5 gal	> 24 hrs	> 24 hrs	> 24 hrs	21 hrs	7.5 hrs	6 hrs
7931	16 oz/5 gal	> 24 hrs	> 24 hrs	12 hrs	23 hrs	9 hrs	7 hrs
7923	8 oz/5 gal	> 24 hrs	12 hrs	8 hrs	20 hrs	10 hrs	8 hrs
7923	16 oz/5 gal	11 hrs	4 hrs	2.5 hrs	13 hrs	7 hrs	5.5 hrs
7925	8 oz/5 gal	5 hrs	2.6 hrs	2.3 hrs	23 hrs	8 hrs	6 hrs
7925	12 oz/5 gal	3 hrs	1 hr	1 hr	20 hrs	6.5 hrs	5 hrs
7925	16 oz/5 gal	3 hrs	0.5 hrs	0.5 hrs	14 hrs	5 hrs	4.5 hrs
7925	20 oz/5 gal	2.5 hrs	0.3 hrs	0.3 hrs	11 hrs	4.5 hrs	4 hrs

7430

Accelerator	Amount	Pot Life @40°F/60%RH	Pot Life @75°F/50%RH	Pot Life @95°F/50%RH	Recoat Time @40°F/60%RH	Recoat Time @75°F/50%RH	Recoat Time @95°F/50%RH
none	none	> 7 days	> 7 days	> 7 days	37 hrs	14 hrs	11 hrs
7931	4 oz/5 gal	> 24 hrs	> 24 hrs	> 24 hrs	24 hrs	12 hrs	10 hrs
7931	8 oz/5 gal	> 24 hrs	> 24 hrs	14 hrs	20 hrs	9 hrs	8 hrs
7931	16 oz/5 gal	> 24 hrs	16 hrs	11 hrs	10 hrs	6 hrs	5 hrs
7923	8 oz/5 gal	> 24 hrs	15 hrs	12 hrs	26 hrs	12 hrs	10 hrs
7923	16 oz/5 gal	20 hrs	10 hrs	8 hrs	20 hrs	13 hrs	9 hrs
7925	8 oz/5 gal	> 24 hrs	> 24 hrs	> 24 hrs	12 hrs	5 hrs	5 hrs
7925	12 oz/5 gal	20 hrs	8 hrs	6.5 hrs	8 hrs	3.5 hrs	3 hrs
7925	16 oz/5 gal	8 hrs	1.7 hrs	1.5 hrs	6 hrs	2.5 hrs	2 hrs

NEOGARD
Div. of JONES-BLAIR

2728 Empire Central - Dallas, Texas 75235 - Phone (214) 353-1600 - Fax (214) 357-7532 - www.neogard.com